

QuattroflowTM EZ-Set

Faster Replacement of Single-Use Pump Chambers

With the EZ-Set, you can change your Single-Use pump chamber in 30 seconds without the use of special tools or torque wrenches. The comfortable handgrip design makes replacing the chamber so easy, even while wearing rubber gloves. EZ-Set reduces downtime between batches, allowing you to spend more time doing what is really important for your work. It's quick, it's easy, and it's safe.

- Easy installation of Single-Use pump chambers
- For Single-Use pump chamber sizes 150 and 1200 (QF150SU / QF1200SU / QF1200SU-CV / QF1200SU-HT)
- Nothing to screw in
- Easy visual guide for proper tightening, no torque wrenches needed
- Available for new pumps by adding "EZ" to the pump code
- Retrofittable on existing standard Quattroflow Single-Use pumps by replacing pump chamber, pressure plate and shaft bearing kit

Order information for a new pump (example):

Pump	Order Number
QF150SU	QF15SU-EZ
QF1200SU	QF12SU3-EZ QF12SU5-EZ

Order information to upgrade from a standard pump to an EZ pump:

Pump	Single-Use Pump Chambers (3 units)	Pressure Plate Including Bolts	Shaft Bearing Kit
QF150SU	QF15DISPP-EZ	PQ15DISKIT-EZ	PSKITWLC155-EZ
QF1200SU	QF12DISPP-EZ	PQ12DISKIT-EZ	PSKITWLC125-EZ

Don't lose more time!

Order your EZ-Set pump or pump upgrade right away!

See how easy and quick EZ-Set really is!



Check out the video at quattroflow.com

Quattroflow EZ-Set

What Was Easy Before, Is Even Easier Now!

Upgrade your QuattroflowTM pump with the EZ-Set pressure plate for Single-Use Quattroflow QF150 and QF1200 sizes



ALMATEC Maschinenbau GmbH
Carl-Friedrich-Gauß-Str. 5
47475 Kamp-Lintfort, Germany
O: +49/2842/961-0, F: +49/2842/961-40
info@almatec.de
quattroflow.com